








Prism Clip-ons For Restricted Head Movement

Prism Applications – Vertically Restricted Neck and Eye Movement



	<p>Base-Up Occupational Prism lowers overhead images into view without the need to arch the neck backward</p>	
	<p>Base-Up Postural Prism lowers images into view. Aids those whose heads cannot be elevated.</p>	
<p>Both Base-Up prism Clip-ons are useful to lower straight and overhead images into view.</p>		
	<p>Base-Down Elevates food and reading material into view if the head or eyes are locked in a forward position and a downward gaze is difficult.</p>	
<p>800.410.1618</p>	<p>Chadwick Optical Inc.</p>	<p>www.chadwickoptical.com</p>